



Temporary Abutment /Esthetic-low Healing Cap

Impression Component

Prosthetic Restoration


Abutment

Titanium Coping

Auxiliary Parts Screw

Mini

Temporary Abutment



H	GH	Type	D	Φ4.0
10	1.0	Hex	371010	
		Non-hex	371020	
	3.0	Hex	371030	
		Non-hex	371040	


\* Tightening torque: 20-30Ncm

Esthetic-low Healing Cap




377040

Manual Impression



L	D	Φ4.0
11	362210	
15	362240	

Digital Impression



L	D	Φ3.8
10	361310	
15	361320	


Manual Impression



Φ4.8
362470
Open tray

Φ4.8
362480
Closed tray

Digital Impression



Φ4.8
362640
Scanbody



361010



361330




362140



362680


Cement Retained Restoration



H	GH	D	Φ4.5
5.5	2.0	Hex	331010
		1.0	331020
		3.0	331030
		4.0	331040
		5.0	331050
7.0	2.0	1.0	331510
		2.0	331520
		3.0	331530
		4.0	331540
		5.0	331550

\* Tightening torque: 30Ncm


Cement Retained Restoration



H	GH	Type	D	Φ4.5
8.0	2.0	A	333010	
		B	333020	
	4.0	A	333040	
		B	333050	


\* Tightening torque: 30Ncm

Screw Retained Restoration



GH	D	Φ4.8
1.0	337010	
2.0	337020	
3.0	337030	
4.0	337040	
5.0	337050	


\* Tightening torque: 30Ncm



A	GH	D	Φ4.8
17°	2.5	Hex	337130
		3.0	337140
		4.0	337150
30°	3.5	3.5	337160
		4.0	337170
		5.0	337180

\* Tightening torque: 25Ncm


Cement Retained Restoration



H	GH	D	Φ4.5	Hex	Φ5.0	Hex	Φ6.0	Hex	Φ7.0
5.5	2.0	3.0	Hex	331060	331160	331310	331410		
			Hex	331070	331170	331320	331420		
			Hex	331080	331180	331330	331430		
			Hex	331090	331190	331340	331440		
			Hex	331100	331200	331350	331450		
7.0	2.0	3.0	Hex	331110	331210	331360	331460		
			Hex	331120	331220	331370	331470		
			Hex	331130	331230	331380	331480		
			Hex	331140	331240	331390	331490		
			Hex	331150	331250	331400	331500		

\* Tightening torque: 35Ncm

Cement Retained Restoration




H	GH	Type	Φ5.0	Hex	Φ6.0
8.0	2.0	A	333070	333130	
		B	333080	333140	
	4.0	A	333100	333160	
		B	333110	333170	

\* Tightening torque: 35Ncm

Regular

Temporary Abutment



H	GH	Type	D	Φ4.5
10	1.0	Hex	371050	
		Non-hex	371060	
	3.0	Hex	371070	
		Non-hex	371080	

\* Tightening torque: 20-35Ncm

Esthetic-low Healing Cap



377040

Manual Impression




L	D	Φ4.0	Φ5.0	Φ6.0
11	362210	362220	362230	
15	362240	362250	362260	

Digital Impression



L	D	Φ4.0	Φ5.0	Φ6.0
11	362300	362320	362340	
14	362310	362330	362350	

Digital Impression



L	D	Φ4.5
10	362570	
15	362580	


Manual Impression



Φ4.8
362470
Open tray

Φ4.8
362480
Closed tray

Digital Impression



Φ4.8
362640
Scanbody




362140




362680

Screw Retained Restoration



GH	D	Φ4.8
1.0	337070	
2.0	337080	
3.0	337090	
4.0	337100	
5.0	337110	

\* Tightening torque: 30Ncm



A	GH	D	Φ4.8
17°	2.5	Hex	337190
		3.0	337200
		4.0	337210
30°	3.5	3.5	337220
		4.0	337230
		5.0	337240

\* Tightening torque: 30Ncm

Overdenture

Healing Cap



325010 325020 325030

Equator abutment number	338010	338050	338250
	338020	338060	338260
	338030	338230	
	338040	338240	
	338210		
	338220		




368010



368020

Retentive Anchor



Gingiva Height	Mini		Regular
	1	2	3
17°	1	338010	338210
	2	338020	338220
	3	338030	338230
30°	4	338040	338240
	5	338050	338250
	6	338060	338260

\* Tightening torque: 35Ncm

SmartBox Male Processing Package 378100



TI Screw

Instrument

Hex Hand Drill



353010 (Short)  
353020 (Long)

Hex Torque Drill



353030 (Short)  
353040 (Long)

Multi Abutment outer Driver



354050

Multi-angled Abutment Carrier




353060

Hex Hand Drill



353010 (Short)  
353020 (Long)

Hex Torque Drill




353030 (Short)  
353040 (Long)

Multi Abutment Outer Driver



354050

Multi-angled Abutment Carrier



353050

Esthetic-low Polishing Protector



357020

Esthetic-low Polishing Protector



357020

Insertion-Extraction Tool



358020